


Hear the colors of sound.

The new Ace hearing instruments. Superb hearing you dont't have to think about.






Ace hearing instruments:

Discreetly sophisticated.


Good things come in small packages. And Ace hearing aids are no exception. Ultra-discreet and packed with features, they're IP67-rated to resist moisture, sweat, dust and dirt for greater reliability.

Ace at a glance:





4



Ace's range of natural colors are designed to match a variety of hair and skin tones.



Enjoy life's true sound.

A noisy bar, a railway station or a street café. Anywhere where there's a high level of ambient noise, communicating with people calls for concentration and effort. Straining to hear voices in difficult situations can, by the end of a long day, be exhausting. It's a situation that even people with normal hearing are familiar with.

The new Ace helps overcome this. Its advanced SpeechMaster function singles out the dominant speaker, reducing extraneous voices and sounds. You hear what you want but not the rest. As a result, hearing becomes effortless and less tiring than normal hearing. All day long. In every situation.





A day out with the kids. And Ace.

Chaos reigns.

At breakfast, the kids are excited about the day. Chatter mingles with the radio while your partner is immersed in the newspaper. Ace brings outstanding performance to the table to make sure you don't miss a thing. Its advanced automatic **SpeechMaster** function singles out the speech while simultaneously engaging **SoundSmoothingTM** to reduce the clattering of cutlery and rustling newspaper.

So you can relax and effortlessly understand every word. With Ace, you concentrate on the conversation, not on hearing.



Catching the breeze.

It's late afternoon and you're sitting at a street café enjoying a cappuccino with friends, the breeze is cool but the discussion at the next table is getting heated. With Ace, it's refreshingly easy. Its automatic **SpeechMaster** function reduces surrounding voices by singling out the discussion at your table while further excluding the heated discussion with directional microphones. It also engages **eWindScreenTM** to tackle the intrusive breeze. You can sit back and enjoy the company.

Echoes of the past.

You're with the children at the museum. Voices are echoing around and mixing with the sound of footsteps on the tiled floors. From behind you, a guide asks if you need any help. This could all be tiring but Ace takes it all in its stride. Its **EchoShield** program all but eliminates sound reflections such as reverberation, delivering a cleaner, clearer sound. The only thing you need to concentrate on is which exhibit to see first.



Made for music lovers.

Ace delivers outstanding quality in every situation - even when you're enjoying music. Whether you're enjoying your favorite music at a concert, listening at home on your own HiFi system, or even on stage performing music yourself, its **HD Music** program gives you fantastic sound quality and a superb music experience. Simply and conveniently activated via the push button or touchControl App, ensuring you enjoy every note.

Tinnitus? Relief is at hand.

With its static presets and ocean wave **tinnitus therapy signals**, Ace provides welcome relief for tinnitus sufferers. Ask your Hearing Care Professional which solution is best for you.

Extra touches for greater convenience.

Complete discretion in the palm of your hands.

Compatible with Android and iOS smartphones, the **touchControl App** lets you change programs, adjust the volume, bass and treble with complete discretion. And here comes the best part: Downloading the app from the Google Play Store or Apple App Store is completely free of charge.



The touchControl App lets you discreetly control Ace via your smartphone.





The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc.

Legal Manufacturer

Signia GmbH Henri-Dunant-Strasse 100 91058 Erlangen Germany

Signia GmbH is a Trademark Licensee of Siemens AG.

Order No. 02704-99C1-7600 | Printed in Germany © 02.2016 Signia GmbH

signia-hearing.com